

Waterproof handheld measuring device for pH / Redox with external probe



**WATER-PROOF
DEVICE AND PLUG CONNECTIONS**

- Serial Interface
- analog output (GMH 5550)
- Data logger- and alarm function (GMH 5550)
- GLP-features (Good Laboratory Practice)
- Robust silicone protection cover
- Big dual display with background illumination
- High resolution (0.001pH / 0.1 mV)
- Incl. calibration protocol

GMH 5530

Waterproof handheld measuring device without electrode

GMH 5550

Waterproof handheld measuring device with analog output and data logger, without electrode

Application:

- Waters measuring, fishkeeping, aquafarming
- Drinking water monitoring, process control, soil measuring
- Food production and monitoring
- Laboratory: Medicine, pharmaceuticals, chemistry
- Quality management

Specification:

Measuring ranges:

pH:	-2.000 ... 16.000 pH
Redox / mV:	-2000.0 ... 2000.0 mV (for hydrogen system DIN38404: -1792 ... +2207 mV _H)
Temperature:	-5.0 ... +150.0 °C 23.0 ... 302.0 °F
rH:	0.0 ... 70.0 rH

Accuracy:

pH:	±0.005 pH
Redox / mV:	±0.05 % FS (mV or mV _H)
Temperature:	±0.2 °C (in the range of -5,0 ... 100,0 °C)
rH:	±0.1 rH



Connections:

pH, Redox:	BNC-female connector, compatible to standard BNC-plugs and waterproof BNC-plugs, additional banana-jack (4 mm) for separate reference electrode, input resistance: 10 ¹² Ohm
Temperature:	2 banana-jacks (4 mm) for temperature probes (Pt1000 or NTC 10K)
Interface / Supply:	4-pole bayonet connector for serial interface and supply (with accessory USB 5100)

Working conditions: -25 ... 50 °C; 0 ... 95 % RH (non condensing)

Display: two 4½ - digit seven-segment display (15 mm and 12 mm)

pH-Calibration:

Automatically:	1-, 2- or 3- point calibration, GREISINGER-Standard-Buffer or Puffer to DIN19266 (A,C,D,F,G)
Manually:	1-, 2- or 3- point calibration

Housing: impact resistant ABS housing with pop-up clip

Protection class: IP65 / IP67

Dimensions: 160 x 86 x 37 mm (H x W x D) incl. protection cover

Weight: 250 g incl. battery and protection cover

Power supply: 2 x AAA-battery (incl. in scope of supply)
power consumption: <1.0 mA

Battery life time: 1000 hours

Functions:

Min / Max Value Memory: highest and lowest measured value is saved

Hold: displayed value gets frozen by keypress

Auto-Hold: automatic freezing of a constant measuring value

Auto Power Off: device is automatically switched off after a selectable period if unused (0 to 120 min, or deactivated)

Additional Display for pH-Electrode and Battery: Bar graph display

Low Battery Display „BAT“

Background illumination: duration adjustable (off, 5 s ... 2 min)

Automatic Temperature Compensation:

There is an automatic temperature compensation (ATC) in the range of 0-105 °C for operation mode "pH" and if a temperature probe is connected. Without connected probe the temperature can be input manually.

pH-Calibration:

The used buffer is detected automatically. The temperature dependency of the buffer is automatically compensated.
Permissible electrodes' data: Asymmetry: ±55 mV / Slope: 45 ... 62 mV/pH
The condition of pH-Electrode is checked at each calibration.
1-, 2- or 3- point calibration with characteristics bend for GREISINGER-Standard-Buffer, buffer to DIN 19266 or manual buffer input

Redox-Measurement (ORP): 2 choices:

"mV" Standard-redox-, ORP or mV- measurement
„mV_H“ Conversion to hydrogen systems according to DIN38404 Teil 6

rH-Measurement: The rH-value is calculated from a measured Redox-value and a manually input pH-value.

	GMH 5530	GMH 5550
Adjustable calibration intervals (GLP)	x	x
Calibration memory (GLP)	-	x (last 16 Calibrations)
Real-time clock	-	x
Analog output	-	0 - 1 V, freely adjustable, connection with 4-pole bayonet connector, Resolution 13 bit, accuracy 0.05 % at nominal temperature
Data logger	-	With measuring point input Recording interval: 1 s ... 1 h Recording period: 416 days at interval 1 h Value memory: cyclic: 10000 data sets; manual: 1000 data sets
Min-/max-alarm	-	Permanent monitoring of alarm boundaries (pH / mV and temperature) 3 alarm conditions - off: Alarm function inactive - on: Alarm report via display, integrated buzzer and interface - no Sound: Alarm report only via display and interface

pH / Redox accessories

Accessories and spare parts:

GMH 55 ES

Supplementary set, including pH-electrode (GE 100 BNC), temperature probe (GTF 55 B), case (GKK 3500), working and calibration set (GAK 1400)

GE 125-BNC


waterproof pH-electrode with integrated Pt1000 temperature sensor incl. waterproof BNC-plug and two banana plugs (p.r.t. page 40)

**GE 117-BNC**

pH-electrode with integrated Pt1000 temperature sensor (p.r.t. page 40)

**GTF 55 B**

Pt1000 temperature immersion sensor for liquids
1 m PVC-cable with two banana plugs

**GE 100-BNC**

pH-electrode (p.r.t. page 40)

**GE 105-BNC**

Redox-electrode (p.r.t. p. 40)

**GAK 1400**

Working and calibration set

PHL 4

ready-to-use buffer solution (pH 4.01 / 25 °C) 250 ml

PHL 7

ready-to-use buffer solution (pH 7.00 / 25 °C) 250 ml

PHL 10

ready-to-use buffer solution (pH 10.01 / 25 °C) 250 ml

**KCL 3 M**

3 mol KCl electrolyte for refill or storage (filled in the protective cap) of electrodes with 3 mol KCl electrolyte. 100 ml plastic vial.

CaCl

1000 ml, solution for measuring the pH value of soil

GRL 100

Pepsin cleansing solution, 100 ml

GRP 100

ORP buffer solution (220 mV at 25 °C), 100 ml

EBS 20M

Software for long-term monitoring (p.r.t. page 66)

GSOFT 3050

Software for operation of logger devices (p.r.t. p. 66)

USB 5100

Electrically isolated interface converter with supply of device via USB

GNG 5 / 5000

Plug-in power supply 5 V DC, suitable for GMH 5000 - series (p.r.t. page 65)

GKK 3500

Device case (394 x 294 x 106 mm) with eggcrate foam and cut-outs for 1 device (p.r.t. page 64)

digital pH-Meter

**GPH 014**

Device ready for use incl. pH-electrode type GE 014 and battery.
(no buffer solutions)

Specification:

Measuring range (device):	0.00 ... 14.00 pH
Resolution:	0.01 pH
Accuracy (device):	0.02 pH ± 1 digit (at nominal temperature = 25 °C)
Input resistance:	10 ¹² Ohm
pH-electrode:	combined measuring and reference electrode type GE 014 with refillable 3 mol-KCl electrolyte, 2-12 pH, 0 ... 60 °C
Calibration:	3 turning knobs for: - temperature compensation 0 ... 90 °C - pH7 value - pH x-value (e.g. 4.0, 10.0, 12.0)
Working temperature:	0 ... 45 °C
Display:	3½-digit LCD display, 13 mm high
Power supply:	9 V battery type JEC 6F22 (incl.)
Battery service life:	approx. 200 operating h
Low battery warning:	automatic; "BAT" displayed in case of low voltage
Dimensions:	approx. 106 x 67 x 30 mm (H x W x D) Impact resistant ABS housing
Weight:	approx. 200 g (incl. battery and electrode)

**GAK 1400**

Working and calibration set

General:**Working and calibration set consisting of:**

5 buffer capsules each GPH4.0, GPH7.0 and GPH10.0, 3 x 100ml-plastic bottle GPF100, 1 x 3 mol KCL-electrolyte KCL3M and 1 x Pepsin-cleaning agent GRL100.

GAK1400 is required if no buffer solutions are existing.

Accessories and spare parts:**GE 014-Cinch**

Spare electrode

GPH 014 GL

Loose device (without accessories)

GE 100-Cinch

Better electrode (0-14 pH, 0-80 °C)

GE 101-Cinch

Injection electrode (2-11 pH, 0-60 °C)

GE 104-Cinch

pH-electrode for low-ion water (as of 25 µS/cm)

GKK 252

Case (235 x 185 x 48 mm) with foam lining

GKK 1100

Case (340 x 275 x 83 mm) with foam lining

GB 9 V

Spare battery

additional accessories p.r.t. p. S. 41