

Application areas: Electrodes

Application	Electrode:									
	GE 100	GE 101	GE 104	GE 108	GE 120	GE 151	GE 170	GE 173	GE 117	GE 125
Sewage								•		
Aquarium water	•			•		•			•	
Emulsions		•	•							
On-site measurements				•					•	•
Fish farming	•		•	•		•			•	•
Galvanic baths								•		
Beverages						•		•		•
Low-ion media			•							
Cosmetics			•							
Food sample		•			•					
Sea water	•	•	•	•	•	•	•	•	•	•
Process chemistry							•	•		
Swimming pool water	•			•				•		
Suspensions		•	•							
Soil testings		•								
Drinking water	•		•					•		
Water-based lacquers			•					•		